

The **Dyna-Form Jupiter** mattress can be used with patients considered to be up to that of "**High Risk**" of pressure ulcer development. Using modern foam designs, castellated on one side to aid pressure ulcer management, whilst providing excellent patient comfort; this mattress can be used without a turning programme, aiding moving & handling procedures. **Jupiter** can work comfortably with any electric profiling bed frame. Its durable design uses only the very best CME foam & fabrics to provide a modern mattress complete with infection control enhancing cover.

Also Available in Pro-Formance™ Strikethrough Resistant Cover.

	Features	Benefits	Technical Specifications	
(1a) (1b)	Castellated foam (with tear duct cut)	1a smaller castellated cells for heels and 1b larger cells sacrally	Product Code Risk Category	MAT/JUP/198/88/15 High Risk
2	Single sided cut surface	Works without a turning programme	Weight Limit Warranty	28 Stone / 180kg 3 Years
3	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear & friction forces, whilst assisting infection control policies	Dimensions	(Length) 198cm x (Width) 88cm x (Height) 15cm — However Variable Sizes Available Upon Request
4	Fully sewn or high frequency welded cover. (NB: welded cover option carries a small price premium)	Prevents fluid ingress & aids decontamination processes (can be laundered up to 95 °C)	Foam Density	34–36kg/m³ Nominal Hardness – 125 Newtons
5	Fully replaceable parts	Ensures great value for money	Fatigue Class Fire Retardancy	Very Severe BS7177; 2008 (Crib 5 Composite)
			Weight	10kg
			Colour	Light Blue
			Hypochlorite	1,000 ppm diluted or 0.1%

Solution

compatible

(or to be used in conjunction with local

infection control guidelines)