

Reynard

Surface Disinfectant Wipes



Surface Disinfecting Wipes

Reynard's surface disinfecting products are alcohol, chlorine and fragrance free to avoid any tainting problems should food be placed in contact with a wiped surface. They contain a very low foaming, biodegradable non-ionic surfactant, which exhibits good detergency for the emulsification of fatty/greasy soils. These surface disinfection products contain a range of Active Biocides designed to target a broad spectrum of micro-organisms.



Directions: Use the required number of wipes.
Dispose of used wipes properly; do not flush.



1 step cleaning & disinfection

Hygienic, easy to use, convenient & efficient. Reynard Premier Detergent & Disinfectant wipes will clean and disinfect most hard surfaces and equipment in one easy step. They are safe to use on mattress covers, beds, tables, work tops, trays, trolleys and other hospital equipment.

Fast and effective

Reynard Premier Detergent & Disinfectant Wipes will kill 99.9% of common germs and have been tested to be effective against a wide range of bacteria. Disposable and soap free.

Cost-effective

Reynard Premier Detergent & Disinfectant wipes will clean and disinfect in one. There is no need for extra products. Reynard wipes eliminates the need for rags, spray bottles and chemicals, thus saving on laundering, chemical and storage costs. The non tainting, low foaming formula does not need to be rinsed from surfaces and the large wipe size of 33 x 22cm ensures you use less wipes for all applications.

Technical Performance:

Bactericidal—MRSA, VRE, E Coli, Pseudomonas Aeruginosa, Staphylococcus aureus, Listeria, Salmonella & Enterococcus hirae,
Mycobactericidal—Mycobacterium avium, Mycobacterium terrae,
Fungicidal—Aspergillus Niger,
Virucidal—Influenza A (H1N1), HIV, Hep B, Vaccina, Norovirus,
Sporicidal—C Diff

These Multi Surface Disinfection products have been tested as per EU Standards, EN1276, EN13727, EN14348, EN13704, EN1650 and EN14476

Product specifications and availability may vary according to country.
Always read the label and product information before use.

Supporting documentation: Full technical specifications are available upon request.
Bespoke solutions: Should your market require more or fewer wipes, lightly scented alternatives or even biodegradable wipes, for example. Let us know. We are more than happy to develop bespoke solutions.

Surface Disinfecting Wipes



One step for cleaning and disinfecting. Reynard Premier Detergent & Disinfectant Wipes are fragrance free to avoid any risk of tainting. They contain a very low foaming, biodegradable non-ionic surfactant, which exhibits good detergency for the emulsification of fatty/greasy soils. These surface disinfection products contain a range of Active Biocides designed to target a broad spectrum of micro-organisms. They are disposable, soap free and effective against 99.9% of common germs.



RHS 216 Premier Detergent & Disinfectant Wipes

Technical specifications

- 33 x 20cm, 40 gsm non-woven
- 50 wipes per pack, 12 packs per carton



RHS 218 Premier Detergent & Disinfectant Wipes

Technical specifications

- 33 x 20cm, 40 gsm non-woven
- 100 wipes per pack, 6 packs per carton



RHS 214 Premier Detergent & Disinfectant Wipes

Technical specifications

- 28 x 25cm, 30 gsm non-woven
- 120 wipes per pack, 4 buckets per carton

Dispenser

- Perspex wall mounted dispenser



RHS 212 Surface Disinfecting Spray / Pour Bottle

Technical specifications

- 1 Litre spray nozzle or pour lid bottle
- 6 bottles per carton



The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. Regarding material compatibility are only possible in particular cases. Our recommendations are without obligation and do not constitute a warranty. They do not preclude the company's own testing for the intended purposes and processes. In this respect we cannot accept any liability.

This complies with our general conditions of sale and supply.

United Kingdom

Boston House
Grove Technology Park
Downsview Road, Wantage
Oxfordshire OX12 9FF
☎ +44 845 263 7335

Australia

Suite 403, 44 Hampden Road
Artarmon
Sydney, NSW 2064
enquiries @reynardhealth.com
☎ 1300 667 199

New Zealand

17 Napier Road
Havelock North
Hawke's Bay
4130 New Zealand
☎ +64 (0)6 650 0709

International

17 Napier Road
Havelock North
Hawke's Bay
4130 New Zealand
☎ +64 (0)6 650 0709

www.reynardhealth.com